


## NEC LOC

### Cervical collar



Bivalve cervical collar designed for use in first aid and emergency situations. It's identified as a solution for the immobilization of the cervical spine in trauma patients.

#### Specific features

- Quick application due to strap closure
- Access for tracheotomy
- Comfortable soft padding

#### Technical data

Size	Small
Materials	PE/Plastizote
Weight	180 g

Class I Med.Device according to Dir. 93/42/CEE

NEC LOC SIZE 6 (SMALL)-  
TWO-PIECE RIGID CERVICAL  
COLLAR YELLOW



JM10040A

Distributed by Spencer Italia S.r.l.

Rev.2 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)