

NEC LOC

Kit 3 cervical collars



Bivalve cervical collars designed for use in first aid and emergency situations. It's identified as a solution for the immobilization of the cervical spine in trauma patients.

Specific Features

- Quick application due to strap closure
- Access for tracheotomy
- Comfortable soft padding

Technical data

| | |
|-------------------|------------------------------|
| Included sizes | Pediatric Small Medium |
| Materials | PE/Plastizote |
| Bag dimensions | 42x34 cm |
| Weight of the kit | 610 ± 10g |

Class I Med.Device according to Dir. 93/42/CEE

NEC LOC KIT THREE SIZES
(4,6,7) WITH BAG
JM10100A



Distributed by Spencer Italia S.r.l.

Rev.1 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it