

# NEC LOC

## Kit 6 cervical collars



Bivalve cervical collars designed for use in first aid and emergency situations. It's identified as a solution for the immobilization of the cervical spine in trauma patients.

### Specific Features

- Quick application due to strap closure
- Access for tracheotomy
- Comfortable soft padding
- The kit is supplied with a practical carry bag

### Technical data

Included sizes	P2 2-6 years
	Sz.4 (Pediatric)
	Sz.5 (Stout)
	Sz.6 (Small)
	Sz.7 (Medium)
	Sz.8 (Large)
Materials	PE/Plastizote

Class I Med.Device according to Dir. 93/42/CEE

NEC LOC KIT SIX SIZES  
(2,4,5,6,7,8) WITH BAG  
JM10101A



Distributed by Spencer Italia S.r.l.

Rev.0 (31/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)