

# SPECTAR

## Extrication cervical collar




Bivalve cervical collar designed for use in first aid and emergency situations. It's identified as a solution for the immobilization of the cervical spine in trauma patients.

### Specific Characteristics

- Anatomical shape of the occipital and mandibular part
- Quick application due to strap closure
- Access for tracheotomy
- Dedicated opening for carotid pulse control
- Comfortable contoured padding

### Technical data

Size	S
Materials	LDPE, PE
Weight	190 g
Dimension A	295 mm
Dimension B	110 mm



### Accessories

- JM30028A – Sternal support

Class I Med.Device according to Dir. 93/42/CEE

SPECTAR FOAM BIVALVE  
CERVICAL COLLAR SZ. S



JM30021A

CND Classification Y060312

Registration number 1213183

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)