

# SPECTAR

## Extrication cervical collar




Spectar cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

### Specific features

- Anatomical shape of the occipital and mandibular part
- Quick application due to strap closure
- Access for tracheotomy
- Dedicated opening for carotid pulse control
- Comfortable contoured padding

### Technical data

Size	S
Materials	LDPE, PE
Weight	190 g
Dimension A	295 mm
Dimension B	110 mm



### Optional accessories

JM30028A – Sternal support

Class I MD compliant with UE Reg. 2017/745

SPECTAR FOAM - CERVICAL  
COLLAR BIVALVE SIZE S

JM30021A

CND Classification Y060312

Registration number 1213183



Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)