

SPECTAR

Extrication cervical collar




Bivalve cervical collar designed for use in first aid and emergency situations. It's identified as a solution for the immobilization of the cervical spine in trauma patients.

Specific Characteristics

- Anatomical shape of the occipital and mandibular part
- Quick application due to strap closure
- Access for tracheotomy
- Dedicated opening for carotid pulse control
- Comfortable contoured padding

Technical data

Size	M
Materials	LDPE, PE
Weight	195 g
Dimension A	295 mm
Dimension B	115 mm



Accessories

JM30028A – Sternal support

Class I Med.Device according to Dir. 93/42/CEE

SPECTAR FOAM – BIVALVE
CERVICAL COLLAR SZ. M

JM30022A

CND Classification



Y060312

Registration number

1213184

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it