

SPECTAR

Extrication cervical collar




Spectar cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

Specific features

- Anatomical shape of the occipital and mandibular part
- Quick application due to strap closure
- Access for tracheotomy
- Dedicated opening for carotid pulse control
- Comfortable contoured padding

Technical data

Size	L
Materials	LDPE, PE
Weight	205 g
Dimension A	335 mm
DimensioneB	125 mm



Optional accessories

JM30028A – Sternal support

Class I MD compliant with UE Reg. 2017/745

SPECTAR FOAM - CERVICAL
COLLAR BIVALVE SIZE L

JM30023A

CND Classification



Y060312

Registration number

1213185

Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it