

SPECTAR

Extrication cervical collar



Spectar cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.


Specific features

- Anatomical shape of the occipital and mandibular part
- Quick application due to strap closure
- Access for tracheotomy
- Dedicated opening for carotid pulse control
- Comfortable contoured padding
- Color coded for quick identification of the size

Technical data

Size	Stout
Materials	LDPE, PE
Weight	185 g
Dimension A	95 mm
Dimension B	55 mm
Dimension C	125 mm
Width front extended valve	355 mm
Width rear extended valve	355 mm

Note: Dimensions subject to tolerances of ± 10 mm



Optional accessories

- JM30028A - Sternal support

Class I MD compliant with UE Reg. 2017/745

SPECTAR FOAM - 2 PIECE CERVICAL
COLLAR SIZE STOUT

JM30024A

CND Classification



Y060312

Registration number

1229067

Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it