

# SPECTAR

## Extrication cervical collar



Spectar cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.


### Specific features

- Anatomical shape of the occipital and mandibular part
- Quick application thanks to the strap closure
- Access for tracheotomy
- Dedicated opening for carotid pulse control
- Comfortable contoured padding
- Color coded for to make size identification easier

### Technical data

|                          |          |
|--------------------------|----------|
| Size                     | XXS      |
| Materials                | LDPE, PE |
| Weight                   | 120 g    |
| Dimension A              | 110 mm   |
| Dimension B              | 55 mm    |
| Dimension C              | 135 mm   |
| Front piece opened width | 245 mm   |
| Rear piece opened width  | 235 mm   |

Note: All dimensions are subject to tolerances of  $\pm 10$  mm



### Optional accessories

JM30028A - Sternal support

Class I MD compliant with UE Reg. 2017/745

SPECTAR FOAM - CERVICAL  
COLLAR BIVALVE SIZE XXS

JM30025A

CND Classification Y060312

Registration number 1281138



Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)