

SPECTAR

Extrication cervical collar



Bivalve cervical collar designed for use in first aid and emergency situations. It's identified as a solution for the immobilization of the cervical spine in trauma patients.

Specific features

- Anatomical shape of the occipital and mandibular part
- Quick application thanks to the strap closure
- Access for tracheotomy
- Dedicated opening for carotid pulse control
- Comfortable contoured padding
- Color coded for to make size identification easier

Technical data

Size	XXS
Materials	LDPE, PE
Weight	120 g
Dimension A	110 mm
Dimension B	55 mm
Dimension C	135 mm
Front piece opened width	245 mm
Rear piece opened width	235 mm

Note: All dimensions are subject to tolerances of ± 10 mm



Optional accessories

JM30028A - Sternal support

Class I Med.Device according to Dir. 93/42/CEE

SPECTAR FOAM - CERVICAL
COLLAR BIVALVE SIZE XXS
JM30025A



CND Classification Y060312
Registration number 1281138

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it