

## SPECTAR

### Kit 6 extrication cervical collars



Spectar cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

#### Specific features

Anatomical shape of the occipital and mandibular part  
 Quick application thanks to the strap closure  
 Access for tracheotomy  
 Dedicated opening for carotid pulse control  
 Comfortable contoured padding  
 Carrying bag in nylon with zip closure  
 The bag has three compartments for a quick identification of the size and is equipped with a carry handle and a hook that allows to hung the bag

#### Technical data

Included sizes	XXS, XS, S, M, L, STOUT
Materials	LDPE, PE
Kit weight	1,25 kg
Bag dimensions	700 x 380 mm
Bag color	Black

#### Optional accessories

JM30028A – Sternal support

Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)

**JM30026A**

SPECTAR FOAM - KIT 6 SIZES W/BAG

SET consisting of the following MD compliant with UE Reg. 2017/745

JM30020A - SPECTAR FOAM - CERVICAL COLLAR BIVALVE SIZE XS	Reg. n. 1213181	CND Y060312
JM30021A - SPECTAR FOAM - CERVICAL COLLAR BIVALVE SIZE S	Reg. n. 1213183	CND Y060312
JM30022A - SPECTAR FOAM - CERVICAL COLLAR BIVALVE SIZE M	Reg. n. 1213184	CND Y060312
JM30023A - SPECTAR FOAM - CERVICAL COLLAR BIVALVE SIZE L	Reg. n. 1213185	CND Y060312
JM30024A - SPECTAR FOAM - 2 PIECE CERVICAL COLLAR SIZE STOUT	Reg. n. 1229067	CND Y060312
JM30025A - SPECTAR FOAM - CERVICAL COLLAR BIVALVE SIZE XXS	Reg. n. 1281138	CND Y060312

Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
**Spencer Italia S.r.l.** Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)