

SPECTAR

Kit 6 Extrication cervical collars



Bivalve cervical collars designed for use in first aid and emergency situations. They are identified as a solution for the immobilization of the cervical spine in trauma patients.

Specific Characteristics

Anatomical shape of the occipital and mandibular part

Quick application thanks to the strap closure

Access for tracheotomy

Dedicated opening for carotid pulse control

Comfortable contoured padding

Carrying bag in nylon with zip closure

The bag has three compartments for a quick identification of the size and is equipped with a carry handle and a hook that allows to hung the bag

Technical data

Included sizes	XXS, XS, S, M, L, STOUT
Materials	LDPE, PE
Kit weight	1,25 kg
Bag dimensions	700x380 mm
Bag color	Black

Accessories

JM30028A – Sternal support

Class I Med.Device according to Dir. 93/42/CEE

SPECTAR FOAM - KIT 6 SIZES
W/BAG

JM30026A

CND Classification



Y060312

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it