

## SPECTAR

### Kit 3 Extrication cervical collars



Bivalve cervical collars designed for use in first aid and emergency situations. They are identified as a solution for the immobilization of the cervical spine in trauma patients.

#### Specific Characteristics

- Anatomical shape of the occipital and mandibular part
- Quick application due to strap closure
- Access for tracheotomy
- Dedicated opening for carotid pulse control
- Comfortable contoured padding
- Carrying bag in nylon with zip closure

#### Accessories

JM30028A - Sternal support

#### Technical data

Included sizes	S,M,L
Materials	LDPE, PE
Kit weight	650 g
Bag dimensions	620x310 mm

Class I Med.Device according to Dir. 93/42/CEE

SPECTAR FOAM – COLLARS  
KIT 3 SZ.(S-M-L) W/BAG

JM30029A

CND Classification



Y060312

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)