

SPECTAR

Kit 3 extrication cervical collars



Spectar cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

Specific features

- Anatomical shape of the occipital and mandibular part
- Quick application due to strap closure
- Access for tracheotomy
- Dedicated opening for carotid pulse control
- Comfortable contoured padding
- Carrying bag in nylon with zip closure

Optional accessories

JM30028A - Sternal support

Technical data

| | |
|----------------|--------------|
| Included sizes | S, M, L |
| Materials | LDPE, PE |
| Kit weight | 650 g |
| Bag dimensions | 620 x 310 mm |

JM30029A

SPECTAR FOAM - KIT 3 SIZES (S-M-L) W/BAG

SET consisting of the following MD compliant with UE Reg. 2017/745

| | | |
|--|--------------------|----------------|
| JM30021A - SPECTAR FOAM - CERVICAL COLLAR BIVALVE SIZE S | Reg. n. 1213183 | CND Y060312 |
| JM30022A - SPECTAR FOAM - CERVICAL COLLAR BIVALVE SIZE M | Reg. n. 1213184 | CND Y060312 |
| JM30023A - SPECTAR FOAM - CERVICAL COLLAR BIVALVE SIZE L | Reg. n. 1213185 | CND Y060312 |

Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it