

BLUE SPLINT

Rigid splint with flexible internal structure



Blue Splint splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

Specific features

- Flexible and moldable core
- Strap closure
- Washable up to 40 °C by removing the inner core
- Soft materials to ensure high comfort

Technical data

Splint	Size ⁽¹⁾ (mm)	Weight (g)
Leg	580 x 450 ± 10	360
Materials	Neoprene, Al, Nylon	

⁽¹⁾ Measures taken with expanded splint and expressed as Length x maximum width

Class I MD compliant with UE Reg. 2017/745

BLUE SPLINT – RIGID SPLINT FOR
LEG

JM80003A

CND Classification



V0804

Registration number

103025

Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it