

## BLUE SPLINT

Rigid splint with flexible internal structure



Immobilization splint for traumatized limbs. Used to limit the damage and pain to the patient related to rescue operations. Ideal for forearm immobilization

### Specific features

- Flexible and moldable core
- Strap closure
- Washable up to 40 ° C by removing the inner core
- Soft materials to ensure high comfort

### Technical data

Splint	Size <sup>(1)</sup> (mm)	Weight (g)
Arm	390x300	206
Materials	Neoprene, Al, Nylon	

<sup>(1)</sup> Measures taken with expanded splint and expressed as Length x maximum width

Class I Med.Device according to Dir. 93/42/CEE

BLUE SPLINT FOREARM

JM80005A



CND Classification

V0804

Rev.1 (27/11/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)