

BLUE SPLINT PRO

Rigid splint with flexible internal structure



Immobilization splint for traumatized limbs. Used to limit the damage and pain to the patient correlated to rescue operations. Ideal for leg immobilization

Specific Characteristics

- Flexible and moldable core
- Quick release buckles for a quick application
- Washable up to 40 ° C by removing the inner core
- Soft materials to ensure high comfort

Technical data

Dimensions	580x450 mm
Materials	Neoprene, Al, Nylon
Weight	510 g

⁽¹⁾ Measures taken with expanded splint and expressed as Length x maximum width

Class I Med.Device according to Dir. 93/42/CEE

BLUE SPLINT PRO LEG - RIGID
SPLINTS WITH FLEXIBLE INTERNAL
STRUCTURE



JM80030A

CND Classification V0804

Registration number 103067

Rev.0 (30/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it