

BLUE SPLINT PRO

Rigid splints with flexible internal structure



Immobilization system for traumatized limbs. Used to limit the damage and pain to the patient correlated to rescue operations.

Specific Characteristics

- Flexible and moldable core
- Quick release fastener
- Suitable to be used as an armrest during intravenous infusions
- Washable up to 40 ° C by removing the inner core
- Soft materials to ensure high comfort

Technical data

Splint	Size ⁽¹⁾ (mm)	Weight (g)
Leg	580x450	510
Arm	590x280	440
Forearm	390x300	270
Wrist	310x240	230
Elbow/ankle	540x280	390
Materials	Neoprene, Al, Nylon	
Weight of the kit	2 kg	
Bag dimensions	620x310 mm	

Standard equipment

JM80008A	Transport bag for splints
----------	---------------------------

Class I Med.Device according to Dir. 93/42/CEE

BLUE SPLINT PRO SET 5 SIZES W/
BAG

JM80035A

CND Classification



V0804

Registration number

103056

⁽¹⁾Measures taken with expanded splint and expressed as Length x maximum width

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it