

# Urrà Splint

## Semi-rigid splint



Semi-rigid splint for traumatized limbs. Used to limit the damage and pain to the patient correlated to rescue operations.

### Specific features

- Low storage space required
- Soft padding to increase patient comfort
- High moldability
- X-Ray compatible

### Technical data

Dimensions	46x11x0,5 cm
Materials	Al, PP
Weight	82 g
<small>Dimensions of length and width subject do tolerances of <math>\pm 0,5</math>cm</small>	

Class I Med.Device according to Dir. 93/42/CEE

URRA' SPLINT – SEMI RIGID SPLINT  
SIZE M

JM90001A

CND Classification



V0804

Registration number

103539

Rev.2 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)