

Urrà Splint

Semi-rigid splint



Semi-rigid splint for traumatized limbs. Used to limit the damage and pain to the patient correlated to rescue operations.

Specific features

- Low storage space required
- Soft padding to increase patient comfort
- High moldability
- X-Ray compatible

Technical data

Dimensions	22x11x0,5 cm
Materials	Al, PP
Weight	43 g
<small>Dimensions of length and width subject do tolerances of $\pm 0,5$cm</small>	

Class I Med.Device according to Dir. 93/42/CEE

URRA' SPLINT – SEMI RIGID SPLINT
sz. S

JM90002A

CND Classification



V0804

Registration number

103550

Rev.2 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it