

Urrà Splint

Semi-rigid splint



Urrà Splint splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

Specific features

- Low storage space required
- Soft padding to increase patient comfort
- High moldability
- X-Ray compatible

Technical data

Dimensions	22 x 11 x 0,5 cm
Materials	Al, PP
Weight	43 g

Dimensions subject to tolerances of ± 10 mm.

Class I MD compliant with UE Reg. 2017/745

URRÀ SPLINT SIZE S - SEMI-RIGID
SPLINTS

JM90002A

CND Classification



V0804

Registration number

103550

Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it