

Urrà Kit

Immobilization kit



Kit of semi-rigid splints for traumatized limbs to be used to limit the damage and pain to the patient correlated to rescue operations.

Specific Characteristics

The kit includes

- 3 Urrà splints
- Tuff cut scissors
- 1 Spencer ICE PACK
- Self adherent flexible bandage

Technical data

Kit dimensions	260x85x240 mm
Materials	Al, PP
Weight	718 g
Dimensions of length and width subject do tolerances of ± 10 mm	

Class I Med.Device according to Dir. 93/42/CEE

URRÀ KIT - IMMOBILIZATION KIT
JM90003A



CND Classification	V0804
Registration number	195296

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it