

# Urrà Splint

## Semi-rigid splint



Semi-rigid splint for traumatized limbs. Used to limit the damage and pain to the patient correlated to rescue operations.

### Specific Characteristics

- Low storage space required
- Soft padding to increase patient comfort
- High moldability

### Technical data

Dimensions	905x110x5 mm
Materials	Al, PP
Weight	145 g

Class I Med.Device according to Dir. 93/42/CEE

URRA' SPLINT CM91X11 ORANGE  
WITHOUT SPENCER LOGO

**JM90005A**



**CND Classification** V0804

**Registration number** 103563

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)