

## MATRIX FIRST AID

Semi rigid case



First aid semi rigid case with reflective insert and transport handles.

### Specific features

- The semi rigid case offers protection to the content
- Reflective insert to increase visibility
- Zip closure
- Central separator with mesh pockets
- Mesh pockets on both sides of the case and adjustable separators for a better organization of the equipment

### Technical data

Case dimensions	450x340x105 ± 10mm
Weight	1,27 ± 0,1 kg
Material	EVA

MATRIX FIRST AID  
MX10301C

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)