

PARAMEDIC

First aid kit



First aid kit for professional users.

Specify characteristics and content	Q.ty
RES-Q-LIGHT PEN	1
UNIVERSAL BLACK - TUFF CUT SCISSORS	1
TWK033 ALUM. TORCH	1
DISPOSABLE NITRIL GLOVES	1
NON STERILE TWEEZERS	1
GMX747 TRIANGULAR BANDAGE	1
S-G GAUZE COMPRESSES STERILE	1
TNT STERILE SHEET	1
G-BANDAGE 5cm X 5 m	1
G-BANDAGE 7cm X 5 m	1
G-BANDAGE 10cm X 5 m	1
GMX510 IDROPHILIC COTTON	1
S-PLASTER SILK BANDAID 2,5 cm x 9,14 m	1
S-PLASTER SILK BANDAID 5 cm x 9,14 m	1
SELF ADH.FLEXIBLE BAND. WHITE h 10cm x 6,3	1
SSP - STRIP BANDAGES 20 PCS DIFFERENT SIZES	1
HYDROGEN PEROXIDE 10VOL 100ML	1
VOMIT BAG	1
T.I.V. SPENCER AUTOMATIC TOURNIQUETTE	1
SPENCER ICE PACK PE BOX 25 PCS	1
URRÀ SPLINT SIZE M - SEMI-RIGID SPLINTS	1
HT ADULT - ISOTHERMAL SHEET	1
FLEX MASK SET WITH ORANGE PVC CASE	1

Technical data	
Dimension bag	365x110x245 ± 20mm
Material bag	Nylon

Med. Dev. CE Marked according to dir. 93/42/CE
PARAMEDIC KIT - PROFESSIONAL FIRST AID KIT

MX10307C

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it