

Haas 30

Motor driven vacuum pump



Suction pump to be installed inside of a suitable pneumatic circuit.

Specific Characteristics

Low weight and dimensions

3 holes for installation

Stripped connection cables:

Red wire – Positive

Blue wire – Negative

The pump is equipped with a safety system against accidental reverse of polarity that allows to maintain flows directionality also if the connection cables are switched

Technical data

Dimensions	145x145x75 mm
Weight	1050 g
Max flow	24 ± 3 l/min
Max vacuum	84 kPa (630mmHg)
Power supply	12V DC
power	72W
Absorption	7A

HAAS 30 - ELECTRIC VACUUM
PUMP 27 LT.MIN. 12V
OX02300B

Distributed by Spencer Italia S.r.l.

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it