

POP

Oxygen distribution line AFNOR



System for wall recessed mounting, with two medical oxygen output units, to be installed inside a medical vehicle on dedicated surfaces.

Specific features

- Simple and safe use guaranteed by the adoption of gas-specific connections and quick release mechanism
- Smooth surfaces for an easy cleaning
- The locking mechanism, allows easy installation

Technical data

Overall dimensions	200x100x86 ±5 mm
Dimensions of the non-recessed part	200x100x36,5 ±5 mm
Connection standard	AFNOR
Materials	Nylon, Al, Brass
Input connection	M ¼" G
Tube length	1,5 m
Operating temperature	-20°C/+60°C
Storage temperature	-15°C/+40°C

Compliant with

UNI EN ISO 9170
NF S 90-116
EN 1789

Ass. Med. Dev. According to Dir. 93/42/CEE

MOD.D002 MEDICAL GAS DISTRIB.

LINE AFNOR 2 WAYS

OX04357C



Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it