

## POP

### Oxygen distribution line AFNOR



System for wall recessed mounting, with two medical oxygen output units, to be installed inside a medical vehicle on dedicated surfaces.

#### Specific features

Simple and safe use guaranteed by the adoption of gas-specific connections and quick release mechanism  
Smooth surfaces for an easy cleaning  
The locking mechanism, allows easy installation

#### Technical data

Overall dimensions	300x100x86 ±5 mm
Dimensions of the non-recessed part	300x100x36,5 ±5 mm
Connection standard	AFNOR
Materials	Nylon, Al, Brass
Input connection	M ¼" G
Operating temperature	-20°C/+60°C
Storage temperature	-15°C/+40°C

#### Compliant with

UNI EN ISO 9170  
NF S 90-116  
EN 1789

Class IIb Med. Dev. According to Dir. 93/42/CEE

POP MODULE 3 OXYGEN SOCKETS

AFNOR

OX04412A



Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)