


Oxygen probe



Oxygen probe for the connection of devices to a gas- specific socket compliant with the standard DIN, installed in systems for the distribution of medical gases.

Specific features

Allows quick connection and disconnection of devices from the sockets
Color coded for easy identification

Technical data

Dimensions	Ø41x72 mm
Connection	G. ¼" F
Material	Nickel-plated brass

Compliant with

DIN 13260-2

Class IIb Med.Device according to Dir. 93/42/CEE

OXYGEN INLET DIN 13260-2

F1/4" GAS

OX05022B

CND Classification



Z120309

Distributed by Spencer Italia s.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it