

Oxygen probe



Oxygen probe for the connection of devices to a gas- specific socket compliant with the standard BS 5682, installed in systems for the distribution of medical gases.

Specific features

Allows quick connection and disconnection of devices from the sockets

Technical data

Thread	F ¼" G
Connection standard	BS 5682
Material	Nickel-plated brass

Compliant with

BS 5682

Class IIb Med.Device according to Dir. 93/42/CEE

OXYGEN INLET BS 5682

F1/4"GAS

OX05026B



Rev.0 (23/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it