

MONITOX

VISUALIZER FOR OXYGEN TANK CONTENT



System for visualization of the residual pressure inside oxygen cylinders to be mounted in fixed position inside the sanitary vehicle.

Specific features

- Easy to install
- Led lights to indicate the residual pressure
- Power on/off button
- Two manometer to be installed on pressure reducers of oxygen tanks

Technical data

Dimensions	85x72x25 ±2 mm
Materials	PE,Nylon, PC
Manometers connection	G ¼"
Power supply	12V cc
Manometers cable length	2 m
Power connection cable	1 m
Power consumption	0,5 A
Weight	420 g

Class I Med. Dev. according to dir. 93/42/CEE

L.E.D. VISUALIZER FOR OXYGEN
TANK CONTENTS MONITOX

OX07110C



Rev.2 (15/01/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it