

## REK PRO

Clamp with aluminum connection for O<sub>2</sub>



Clamp to be used for connection and mounting of devices for oxygen subministration.

### Specific features

- Extremely compact and robust
- Quick mounting system
- High resistance to involuntary sideways movements thanks to the use of inserts with high adhesion at the points of contact with the rail bar
- Can be used with rail bars complying with the standard UNI EN 19054

### Technical data

Dimensions	67x67x90 mm
Materials	Nylon, Al
Weight	200 ± 5 g

Class I Med.Device according to Dir. 93/42/CEE

REK PRO CLIP WITH  
ALLUMINIUM SUPPORT FOR O<sub>2</sub>  
**OX07116A**



CND Classification	V9099
Registration number	455156

Rev.0 (25/07/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)