

## REK-OX

Manual switcher



Manual exchanger intended for use with Spencer oxygen delivery systems. The device is designed for use with Rek-Ox systems and allows manual selection of the supply oxygen tank simply by turning the knob.

### Specific features

Maximum integration with the distribution system  
 The selection is intuitive thanks to the big sized identification numbers and to the colour coded positions  
 Compact design achieved thanks to a very small pneumatic valve  
 High structural strength thanks to the one-piece aluminum body

### Technical data

Dimensions	90x70x65 mm
Materials	Al, PP, ABS
Input and output connections	M 1/8
Working pressure	400 <sub>0</sub> <sup>+100</sup> kPa
Weight	0,5 kg

Class IIa Med. Dev. According to Dir. 93/42/CEE

REK-OX SWITCH  
 BOTTLE.1/BOTTLE.2  
**OX07202**



CND Classification	V9099
Registration number	1271659

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)