

# Rek OX

## Oxygen distribution line



Indicative image only, it can differ from the actual product



System for wall mounting, with two medical oxygen output units, to be installed inside a medical vehicle on dedicated surfaces.

### Specific features

Simple and safe use guaranteed by the adoption of gas-specific connections and quick release mechanism  
Smooth surfaces for an easy cleaning

### Technical data

Connection standard	BS
Materials	Nylon, Al, Brass

### Compliant with

BS 5682

Class IIb Med. Dev. According to Dir. 93/42/CEE

REK-OX 2 SOCKETS O2 BS  
OX07918



Rev.0 (23/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)