

OX-VA

Oxygen tank with integrated pressure reducer and flowmeter



Oxygen tank with flowmeter to be used in oxygen therapy.

Specific features

- 2 liter oxygen tank with integrated flowmeter
- Can contain pressurized O₂ up to 200 bar
- High autonomy with reduced dimensions
- Flow adjustment from 0 to 15 l/min

Technical data

Dimensions	ø140x450 ± 5 mm
O ₂ Volume at 200bar	400 l
Connection standard	DIN 477 N°6
Output pressure	0 ÷ 2,8 bar
Output	Hose holder ø6 mm
Flow adjustment	From 0 to 15 l/min
Flow precision	±10% or ± 0,5 l/min
Manometer	0÷315 bar, cl 2,5, ø40 mm
Materials	Steel, PE, Chromed brass, Al, PC

Device marked CE according to dir. 93/42/CE

OX. BOTTLE 2L STANDARD DIN
W/REDUCER+PROTECT

OX10017E

Distributed by Spencer Italia s.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it