

OX-T

Oxygen tank



Storage vessel to contain medical oxygen in gaseous state.

Specific Characteristics

Painted steel body
 Steel disk marked for medical use
 Accompanied by certificate of conformity
 Can contain pressurized O₂ up to 200 bar

Technical data

Diameter	110 mm
Height	390 mm
Liters	2l
O ₂ volume at 200 bar	400 l
Connection	DIN
Weight	3 kg

OXYGEN BOTTLE 2lt DIN W/
 VALVE/RING
 OX10019E



Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it