

Oxygen tank with flowmeter and oxygen socket



Indicative image



Oxygen tank equipped with flowmeter and oxygen socket.

Specific features

5 liters oxygen tank with integrated protection and base to increase stability
Can contain O₂ up to 200 bar

Technical data

O ₂ Volume at 200bar	1000 l
Input	NF (ø22,91x1,814) type A
Socket (output)	AFNOR
Socket output flow	140 ± 10 l/min
Output pressure	3,5 ± 0,5 bar
Flowmeter adjustments	0-1-1.5-2-3-4-6-9-12-15 l/min
Flowmeter output	Hose holder ø7mm
Manometer scale	0-300 bar
Materials	Steel, PE, Chromed brass, Al, PC

Device subject to CE marking according to dir. 93/42/CE

OXYGEN TANK 5lt NF W/FLOWM.

AND NF SOCKET

OX10058B

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it