

OX-T

Oxygen tank



Oxygen tank with protection.

Specific features

5 liters oxygen tank with integrated protection and large base to increase stability
 Can contain pressurized O₂ up to 200 bar
 High autonomy in reduced dimensions

Technical data

Diameter	140 mm
Height	565 mm
Capacity	5l
O ₂ Volume at 200 bar	1000 l
Connection	BS (G5/8" type A)

OXYGEN TANK BS W/VALVE
 AND PROTECTION
OX10063E



Rev.0 (01/12/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it