

OX-Gr

Oxygen tank with pressure reducer
Flowmeter and humidifier



Oxygen tank complete of pressure reducer, flowmeter and humidifier intended for the treatment of respiratory failure.

Specific Characteristics

10 liters oxygen tank with integrated protection and large base to increase stability
Can contain pressurized O₂ up to 200 bar
High autonomy in reduced dimensions
Piston pressure reducer with sintered filter
Humidifier with high transparency to ensure an easy reading of the level

Technical data

Tank dimensions	ø140x940 ± 5 mm
Protection diameter	140 ± 5 mm
Assembly dimensions	970x140x260 ± 10mm
O ₂ volume at 200bar	2000 l
Reducer inlet	NF
Reducer output pressure	3,5/4 bar
Flowmeter regulation	From 0 to 15 l/min
Humidifier indicator of maximum	130 ml
Humidifier output	Hose ø7mm
Materials	Steel, PE, Chromed brass, Al, PC
Overall weight	15,3 ± 0,3 kg

Device subject to CE marking according to the dir.
93/42/CE

MOD.G012 OX-GR TANK 10L
W/RED./FLOW./HUMIDIF.AFNO
R
OX10066D

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it