

OX-T

Oxygen tank



The image can differ from the actual product

Oxygen tank with protection.

Specific features

5 liters oxygen tank with integrated protection and large base to increase stability
 Can contain pressurized O₂ up to 200 bar
 High autonomy in reduced dimensions

Technical data

Tank dimensions	ø140x565±10 mm
O ₂ volume at 200bar	1000 l
Connection	DIN (G 3/4" type B)
Weight	8,3 ± 0,2 kg

OXYGEN TANK 5L DIN WITH
 PROTECTION
OX10068



Distributed by Spencer Italia s.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza(PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it