

Oxygen tank with flowmeter and oxygen socket



Indicative image



Oxygen tank equipped with flowmeter and oxygen socket.

Specific features

10 liters oxygen tank with integrated protection and base to increase stability
Can contain O₂ up to 200 bar

Technical data

O ₂ Volume at 200bar	2000 l
Input	UNI (W 21,7x1/14") Type B
Socket (output)	UNI
Socket output flow	140 ± 10 l/min
Output pressure	3,5 ± 0,5 bar
Flowmeter adjustments	0-1-1.5-2-3-4-6-9-12-15 l/min
Flowmeter output	Hose holder ø7mm
Manometer scale	0-300 bar
Materials	Steel, PE, Chromed brass, Al, PC

Device subject to CE marking according to dir. 93/42/CE

10lt OXYGEN TANK W/FLOWM.
AND SOCKET UNI
OX10098B

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it