

O-Tube

Oxygen tube



Oxygen tube for pneumatic connection of devices included in oxygen distribution lines.

Specific features

Colour and marking clearly identifies the use with oxygen
 Suitable for applications at low pressure (max 14 bar)
 DEHP FREE

Technical data

Internal diameter	6,7 mm
External diameter	12,7 mm
Materials	PVC with fabric reinforcement
Burst pressure	56 bar/23°C

NON-TOXIC TUBE FOR
 MEDICAL OXYGEN Ø12,7x6,7
 EN ISO 5359
OX10104B

Rev.1 (15/01/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it