

OX-T

Oxygen tank



Oxygen tank to contain compressed medical oxygen gas.

Specific features

10 liters oxygen tank with integrated protection and large base to increase stability
 Can contain pressurized O₂ up to 200 bar
 High autonomy in reduced dimensions

Technical data

Diameter	140 mm
Height	940 mm
Liters	10l
O ₂ Volume at 200 bar	2000 l
Connection	BS
Weight	14,4 ± 0,2 kg

O₂ OXYGEN TANK 10lt BS
 OX10112C



Rev.0 (22/09/2020)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it