

OX-T

Oxygen tank



Oxygen tank with protection.

Specific features

10 liters oxygen tank with integrated protection and large base to increase stability
 Can contain pressurized O₂ up to 200 bar
 High autonomy in reduced dimensions

Technical data

Tank dimensions	ø140x940 ± 5 mm
Protection diameter	140 ± 5 mm
O ₂ volume at 200bar	2000 l
Connection	DIN (G 3/4" type B)
Weight	14,4 ± 0,2 kg

BOMBOLA O₂ 10lt DIN TED.
 CON VALVOLA E PROTEZ.
OX10117



Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it