

OXYPACK 2 KIT

Oxygen tank carrier



Tank carrier with Oxygen tank complete of pressure reducer and manual resuscitator with mask.

Specific features

- 2 liters oxygen tank
- Can contain pressurized O₂ up to 200 bar
- High autonomy in reduced dimensions
- Pressure reducer with sintered filter
- Orientable hose to connect oxygen tube easier
- Equipped with disposable resuscitation bag and mask

Technical data

Dimensions	150x290x500 ± 10 mm
Resuscitator volume	1600 ml
Mask material	PVC
O ₂ volume at 200bar	400 l
Connection standard	UNI
Output pressure	3,5/4 bar
Adjustable Oxygen Flow	0- 1-1,5-2-3-4-6-9-12-15 l/min
Pressure reducer output	Hose ø7mm
Overall weight	4,5 ± 0,1 kg

Med.Dev. CE marked according to the dir. 93/42/CE

OXYPACK 2 KIT WITH OXYGEN

TANK 2L AND PRES. REDUC.

OX10203B

Rev.3 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it