

OXYPACK KIT

Oxygen tank with pressure reducer
And oxygen mask



Oxygen tank complete of pressure reducer and oxygen mask.

Specific Characteristics

2 liters oxygen tank
Can contain pressurized O₂ up to 200 bar
High autonomy in reduced dimensions
Pressure reducer with sintered filter

Technical data

Tank dimensions	ø110x390 ± 5 mm
O ₂ volume at 200bar	400 l
Connection standard	DIN
Output pressure	3,5/4 bar
Oxygen Flow	6 l/min
Pressure reducer output	Hose ø7mm

Med.Dev. CE marked according to the dir. 93/42/CE

MOD.D011 OXYPACK W/TANK
2LT DIN/PRESSURE
RED./OXYMASK
OX10204

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it