

## OXYPACK 1 KIT

Oxygen tank with pressure reducer



Oxygen tank complete of pressure reducer.

### Specific Features

- 5 liters oxygen tank
- Can contain pressurized O<sub>2</sub> up to 200 bar
- High autonomy in reduced dimensions
- Easy to carry thanks to the shoulder strap
- Wide pocket with zip closure to store accessories
- Hole with metal ring inside the pocket allows to pass an oxygen tube and to connect it to an oxygen mask
- Pressure reducer with pressure gauge 0-300 bar with sintered filter and integrated flowmeter

### Technical data

Bag Dimensions <sup>(without tank)</sup>	640x270 ± 20mm
Side pocket dimensions	370x125x80 ± 5 mm
Bag weight	600 ± 10 g
Tank dimensions	Ø140x565 mm
Tank weight	8,3 ± 0,1 kg
O <sub>2</sub> volume at 200bar	1000 l
Connection standard	UNI
Output pressure	3,5/4 bar
Adjustable Oxygen Flow	0- 1-1,5-2-3-4-6-9-12-15 l/min
Pressure reducer weight	720 ± 10g
Pressure reducer output	Hose ø7mm

Med.Dev. CE marked according to the dir. 93/42/CE

OXYPACK 1 KIT WITH OXYGEN  
CYLINDER 5 LITERS AND  
PRESSURE REDUCER  
**OX10205B**

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)