

OXYMASK 150

High concentration oxygen masks



Disposable face mask for oxygen therapy.

Specific Characteristics
Adjustable elastic
Semi transparent shell to improve patient monitoring
Anti-crush tube
Flexible aluminum nose clip to increase stability
Comes with reservoir
Latex Free

Technical data	
Material	PVC
Tube length	2 m

Class IIa Med.Device according to Dir. 93/42/CEE	
ADULT OXYGEN MASK HIGH CONCENTRATION WITH RESERVO	
OX20104A	
CND Classification	R03010201
Registration number	104794



Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it