

OXYMASK 100

Variable concentration oxygen mask



Masks for oxygen therapy with venturi system for variation of oxygen concentration.

Specific features
Adjustable elastic
Semi transparent shell to improve patient monitoring
Anti-crush tube
Flexible aluminum nose clip to increase stability
Latex Free

Technical data	
Materials	PVC/PP
Tube length	2 m
Corrugated Tube length	15 cm
Obtainable oxygen concentrations	24/28/31/35/40/50/60%

Class IIa Med.Device according to Dir. 93/42/CEE

PEDIATRIC OXYGEN MASK
VARIABLE CONCENTRATION
OX20106A



Rev.0 (27/07/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it