

# OXYMASK 150

High concentration oxygen masks



Pediatric disposable high concentration face mask for oxygen therapy.

### Specific Characteristics

Adjustable elastic  
Semi transparent shell to improve patient monitoring  
Anti-crush tube  
Flexible aluminum nose clip to increase stability  
Comes with reservoir

**Latex Free**

### Technical data

Material	PVC
Tube length	2 m

Class IIa Med.Device according to Dir. 93/42/CEE

PEDIATRIC OXYGEN MASK HIGH  
CONCENTRATION W/RESERVOIR  
**OX20107A**



<b>CND Classification</b>	R03010201
<b>Registration number</b>	104796

Rev.0 (16/04/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[spencer.it](http://spencer.it)